

SPECIMEN REQUIREMENT FOR 24HOUR URINE TESTING

1. COLLECT THE 24 HOUR URINE IN PROPER SPECIMEN CONTAINER, PROVIDED BY GARCIA LABORATORY
2. USE SMALL ORANGE TOP STERILE URINE CONTAINER FOR SPECIMEN TRANSPORT.
3. FILL ORANGE TOP STERILE URINE CONTAINER WITH A SAMPLE FROM THE 24 HOUR URINE CONTAINER.
4. PROPERLY LABEL ORANGE TOP STERILE URINE CONTAINER WITH PATIENT NAME, INMATE NUMBER AND/OR DOB, AND THE TOTAL VOLUME OF THE 24 HOUR URINE CONTAINER.
5. PLEASE DO NOT SEND 24 HOUR URINE CONTAINER IN THE UPS SHIPMENT. WE ONLY NEED THE TOTAL VOLUME COLLECTED AND THE ORANGE URINE SAMPLE CUP WITH THE SPECIMEN.